




ARRB1 Antibody

Product Code	CSB-PA000888
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49407
Immunogen	Synthesized peptide derived from C-terminal of Human ARRB1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ARRB1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ARR1; arrestin 2; arrestin beta 1; beta-arrestin-1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARRB1

Image



Western blot analysis of extracts from Jurkat cells and HUVEC cells, using ARRB1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.