




CABLES1 Antibody

Product Code	CSB-PA138718
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDN4
Immunogen	Synthesized peptide derived from internal of Human Ik3-1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total Ik3-1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cdk5 and Abl enzyme substrate 1; FLJ35924;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CABLES1

Image



Western blot analysis of extracts from LOVO cells, using Ik3-1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.