





CABLES2 Antibody

Product Code	CSB-PA943938
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTV7
Immunogen	Synthesized peptide derived from internal of Human Ik3-2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total Ik3-2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CDK5 and ABL1 enzyme substrate 2; Interactor with CDK3 2; Ik3-2; CABLES2; C20orf150
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CABLES2

Image



Western blot analysis of extracts from Jurkat cells, using Ik3-2 antibody.



Immunofluorescence analysis of A549 cells, using Ik3-2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.