





CDH17 Antibody

Product Code	CSB-PA212083
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12864
Immunogen	Synthesized peptide derived from internal of Human CDH17.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CDH17 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CAD17; CADH; Intestinal peptide-associated transporter HPT-1; LI-cad; LI-cadherin
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDH17

Image



Western blot analysis of extracts from Jurkat cells and HUVEC cells, using CDH17 antibody.



Immunofluorescence analysis of HeLa cells, using CDH17 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.