




TMEM30B Antibody

Product Code	CSB-PA194152
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3MIR4
Immunogen	Synthesized peptide derived from internal of Human TMEM30B.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TMEM30B protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cell cycle control protein 50B; Transmembrane protein 30B; TMEM30B; CDC50B; CC50B
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TMEM30B

Image



Western blot analysis of extracts from K562 cells, using TMEM30B antibody.