





CDCA4 Antibody

Product Code	CSB-PA597809
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXL8
Immunogen	Synthesized peptide derived from internal of Human CDCA4.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total CDCA4 protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cell division cycle-associated protein 4; Hematopoietic progenitor protein; HEPP;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDCA4

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue using CDCA4 antibody.



Immunofluorescence analysis of HepG2 cells,
using CDCA4 antibody.