




CDCA7 Antibody

Product Code	CSB-PA226976
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWT1
Immunogen	Synthesized peptide derived from internal of Human CDCA7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CDCA7 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cell division cycle-associated protein 7; Protein JPO1; JPO1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDCA7

Image



Western blot analysis of extracts from K562 cells, using CDCA7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.