




CDK5RAP3 Antibody

Product Code	CSB-PA044056
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JB5
Immunogen	Synthesized peptide derived from internal of Human CDK5RAP3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CDK5RAP3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CDK5 regulatory subunit-associated protein 3; CDK5 activator-binding protein C53; Protein HSF-27; IC53; MSTP016
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDK5RAP3

Image



Western blot analysis of extracts from 293 cells, using CDK5RAP3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.