




CNTD2 Antibody

Product Code	CSB-PA081295
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8S5
Immunogen	Synthesized peptide derived from internal of Human CNTD2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CNTD2 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cyclin N-terminal domain-containing protein 2; CNTD2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNTD2

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using CNTD2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.