





COL9A3 Antibody

Product Code	CSB-PA833813
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14050
Immunogen	Synthesized peptide derived from internal of Human Collagen IX α 3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Collagen IX α 3 protein.
Tested Applications	ELISA, IHC, IF; IHC:1:50-1:100, IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Collagen alpha-3(IX) chain ; CO9A3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COL9A3

Image



Immunohistochemistry analysis of paraffin-embedded human tonsil tissue using Collagen IX α 3 antibody.



Immunofluorescence analysis of NIH/3T3 cells, using Collagen IX α 3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.