




COL11A2 Antibody

Product Code	CSB-PA634265
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13942
Immunogen	Synthesized peptide derived from internal of Human Collagen XI α 2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total Collagen XI α 2 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Collagen alpha-2(XI) chain; COBA2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COL11A2

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using Collagen XI α 2 antibody.