





# COL20A1 Antibody

<b>Product Code</b>	CSB-PA236472
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9P218
<b>Immunogen</b>	Synthesized peptide derived from C-terminal of Human Collagen XX $\alpha$ 1.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Specificity</b>	The antibody detects endogenous levels of total Collagen XX $\alpha$ 1 protein.
<b>Tested Applications</b>	ELISA, WB, IHC, IF; WB: 1:500-1:3000, IHC: 1:50-1:100, IF: 1:100-1:500
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Collagen alpha-1(XX) chain [Precursor]; COKA1;
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	COL20A1


## Image



Western blot analysis of extracts from HUVEC cells, using Collagen XX  $\alpha$ 1 antibody.



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue using Collagen XX  $\alpha$ 1 antibody.



Immunofluorescence analysis of NIH/3T3 cells, using Collagen XX  $\alpha$ 1 antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.