





COL20A1 Antibody

Product Code	CSB-PA236472
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P218
Immunogen	Synthesized peptide derived from C-terminal of Human Collagen XX α 1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total Collagen XX α 1 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Collagen alpha-1(XX) chain [Precursor]; COKA1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COL20A1


Image



Western blot analysis of extracts from HUVEC cells, using Collagen XX α 1 antibody.



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue using Collagen XX α 1 antibody.



Immunofluorescence analysis of NIH/3T3 cells, using Collagen XX α 1 antibody.