




COX19 Antibody

Product Code	CSB-PA562512
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q49B96
Immunogen	Synthesized peptide derived from internal of Human COX19.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total COX19 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome c oxidase assembly protein COX19; Hcox19; COX19;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COX19

Image



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using COX19 antibody.