




COX6C Antibody

Product Code	CSB-PA906089
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09669
Immunogen	Synthesized peptide derived from internal of Human COX6C.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total COX6C protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 1.9.3.1; cytochrome c oxidase polypeptide VIc;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COX6C

Image



Western blot analysis of extracts from HUVEC cells, using COX6C antibody.



Immunofluorescence analysis of HepG2 cells,
using COX6C antibody.