





COX7A2P2 Antibody

Product Code	CSB-PA192594
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60397
Immunogen	Synthesized peptide derived from internal of Human COX7S/A2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total COX7S/A2 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome c oxidase subunit VIIa 3; mitochondrial [Precursor]; COX7AP2; COX7A3; COX7AL2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COX7A2P2

Image



Western blot analysis of extracts from rat heart cells, using COX7S/A2 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using COX7S/A2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.