




CYP1A1/CYP1A2 Antibody

Product Code	CSB-PA153187
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04798
Immunogen	Synthesized peptide derived from internal of Human Cytochrome P450 1A1/2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Cytochrome P450 1A1/2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AHH; AHRR; CP11; CYP1; CYPIA1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP1A1/CYP1A2

Image



Western blot analysis of extracts from COS-7 cells, RAW264.7 cells and HT-29 cells, using Cytochrome P450 1A1/2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.