




CYP2A13 Antibody

Product Code	CSB-PA833745
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16696
Immunogen	Synthesized peptide derived from internal of Human Cytochrome P450 2A13.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 2A13 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome P450 2A13; EC 1.14.14.1; CYP11A13; CP2AD;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2A13

Image



Western blot analysis of extracts from HeLa cells, 293 cells and A549 cells, using Cytochrome P450 2A13 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.