





CYP2C8 Antibody

Product Code	CSB-PA694882
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10632
Immunogen	Synthesized peptide derived from internal of Human Cytochrome P450 2C8.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 2C8 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome P450 2C8; EC 1.14.14.1; CYP11C8; P450 form 1; P450 MP-12/MP-20
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2C8

Image



Western blot analysis of extracts from HuvEc cells, using Cytochrome P450 2C8 antibody.



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using Cytochrome P450 2C8 antibody.