




CYP2C18/CYP2C8/CYP2C9/CYP2C19 Antibody

Product Code	CSB-PA834244
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33260
Immunogen	Synthesized peptide derived from internal of Human Cytochrome P450 2C8/9/18/19.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 2C8/9/18/19 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CP2CI; CPCI; CYP1IC18; EC 1.14.14.1; P450-6B/29C
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2C18/CYP2C8/CYP2C9/CYP2C19

Image



Western blot analysis of extracts from LOVO cells, using Cytochrome P450 2C8/9/18/19 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.