




CYP2D6 Antibody

Product Code	CSB-PA263930
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10635
Immunogen	Synthesized peptide derived from internal of Human Cytochrome P450 2D6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 2D6 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome P450 2D6; EC 1.14.14.1; CYP2D6; P450-DB1; Debrisoquine 4-hydroxylase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2D6

Image



Western blot analysis of extracts from HT-29 cells, using Cytochrome P450 2D6 antibody.