





CYP2E1 Antibody

Product Code	CSB-PA961928
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05181
Immunogen	Synthesized peptide derived from C-terminal of Human Cytochrome P450 2E1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 2E1 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CP2E1; CYP2E; CYP11E1; Cytochrome P450 2E1; P450-J
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2E1


Image



Western blot analysis of extracts from K562 cells and HepG2 cells, using Cytochrome P450 2E1 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using Cytochrome P450 2E1 antibody.



Immunofluorescence analysis of HepG2 cells, using Cytochrome P450 2E1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.