




CYP4F2 Antibody

Product Code	CSB-PA916473
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78329
Immunogen	Synthesized peptide derived from N-terminal of Human Cytochrome P450 4F2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 4F2 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome P450 4F2; EC 1.14.13.30; CYP1V2; Leukotriene-B(4) omega-hydroxylase; Leukotriene-B(4) 20-monooxygenase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP4F2

Image



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using Cytochrome P450 4F2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.