





CYP4X1 Antibody

Product Code	CSB-PA071999
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N118
Immunogen	Synthesized peptide derived from internal of Human Cytochrome P450 4X1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 4X1 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CP4X1; CYP4X1; Cytochrome P450 4X1; EC 1.14.14.1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP4X1

Image



Western blot analysis of extracts from COS7 cells, using Cytochrome P450 4X1 antibody.



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using Cytochrome P450 4X1 antibody.