




DNAJB4 Antibody

Product Code	CSB-PA185290
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UDY4
Immunogen	Synthesized peptide derived from internal of Human DNAJB4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total DNAJB4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DJB4; DnaJ homolog subfamily B member 4; DNAJW; DNJB4; Heat shock 40 kDa protein 1 homolog
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DNAJB4

Image



Western blot analysis of extracts from HepG2 cells, using DNAJB4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.