




GBP1 Antibody

Product Code	CSB-PA695454
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32455
Immunogen	Synthesized peptide derived from internal of Human GBP1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total GBP1 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	GBP-1; Guanine nucleotide-binding protein 1; Interferon-gamma inducible protein MAG-1; Interferon-induced guanylate-binding protein 1; MAG-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GBP1

Image



Western blot analysis of extracts from HUVEC cells, using GBP1 antibody.



Immunofluorescence analysis of MCF-7 cells, using GBP1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.