





HSPB2 Antibody

Product Code	CSB-PA052298
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16082
Immunogen	Synthesized peptide derived from internal of Human HSPB2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HSPB2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DMPK-binding protein; heat-shock protein beta-2; HSP27; MKBP;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSPB2

Image



Western blot analysis of extracts from HepG2 cells and HeLa cells, using HSPB2 antibody.



Immunofluorescence analysis of A549 cells, using HSPB2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.