




IL17RA Antibody

Product Code	CSB-PA058027
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96F46
Immunogen	Synthesized peptide derived from internal of Human IL17RA.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total IL17RAprotein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	IL-17 receptor; IL-17RA; IL17RA; interleukin-17 receptor;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL17RA

Image



Western blot analysis of extracts from LOVO cells, using IL17RA antibody.