




IL20RB Antibody

Product Code	CSB-PA125791
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXL0
Immunogen	Synthesized peptide derived from internal of Human IL20RB.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total IL20RB protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DIRS1; FNDC6; I20RB; IL-20R-beta; IL-20R2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL20RB

Image



Western blot analysis of extracts from HeLa cells, using IL20RB antibody.