




MAPK9 Antibody

Product Code	CSB-PA900267
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C12684
Immunogen	Synthesized peptide derived from internal of Human MAPK9.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MAPK9 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 9; EC 2.7.11.24; Stress-activated protein kinase JNK2; c-Jun N-terminal kinase 2; JNK-55
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK9

Image



Western blot analysis of extracts from A549 cells and RAW264.7 cells, using MAPK9 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.