




# PRKAB1 Antibody

<b>Product Code</b>	CSB-PA695332
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y478
<b>Immunogen</b>	Synthesized peptide derived from N-terminal of Human PRKAB1.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous levels of total PRKAB1 protein.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:3000
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	5-AMP-activated protein kinase subunit beta-1; AMPK beta-1 chain; AMPKb; PRKAB1; AMPK
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PRKAB1

## Image



Western blot analysis of extracts from RAW264.7 cells, treated with TNF (20ng/ml, 5mins), using PRKAB1 antibody.