




PRKAB1 Antibody

Product Code	CSB-PA695332
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y478
Immunogen	Synthesized peptide derived from N-terminal of Human PRKAB1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PRKAB1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	5-AMP-activated protein kinase subunit beta-1; AMPK beta-1 chain; AMPKb; PRKAB1; AMPK
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKAB1

Image



Western blot analysis of extracts from RAW264.7 cells, treated with TNF (20ng/ml, 5mins), using PRKAB1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.