




BARD1 Antibody

Product Code	CSB-PA138883
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99728
Immunogen	Synthesized peptide derived from N-terminal of Human BARD1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BARD1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BRCA1-associated RING domain protein 1; BARD-1; BARD1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BARD1

Image



Western blot analysis of extracts from HeLa cells, using BARD1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.