





BAX Antibody

Product Code	CSB-PA552355
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55269
Immunogen	Synthesized peptide derived from internal of Human BAX.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total BAX protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BAX protein; cytoplasmic isoform delta; BAX;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAX

Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using BAX antibody.



Immunofluorescence analysis of A549 cells,
using BAX antibody.