





BCL2L2 Antibody

Product Code	CSB-PA093718
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92843
Immunogen	Synthesized peptide derived from internal of Human BCLW.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total BCLW protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Apoptosis regulator Bcl-W; Bcl-2-like 2 protein; BCL2L2; BCLW; KIAA0271
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BCL2L2

Image



Western blot analysis of extracts from COLO cells, using BCLW antibody.



Immunofluorescence analysis of HepG2 cells, using BCLW antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.