





CIDEB Antibody

Product Code	CSB-PA596246
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHD4
Immunogen	Synthesized peptide derived from internal of Human CIDEB.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total CIDEB protein
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	cell death activator CIDE-B; cell death-inducing DFFA-like effector B;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CIDEB


Image



Western blot analysis of extracts from RAW264.7 cells and HT-29 cells, using CIDEB antibody.



Immunohistochemistry analysis of paraffin-embedded human cervix carcinoma tissue using CIDEB antibody.



Immunofluorescence analysis of COS-7 cells, using CIDEB antibody.