



DCC Antibody

Product Code	CSB-PA081145
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43146
Immunogen	Synthesized peptide derived from internal of Human DCC.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total DCC protein.
Tested Applications	ELISA, IF; IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	colorectal cancer suppressor; deleted in colorectal carcinoma; netrin receptor DCC; tumor suppressor protein DCC;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DCC

Image



Immunofluorescence analysis of LOVO cells, using DCC antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.