




LIG4 Antibody

Product Code	CSB-PA869550
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49917
Immunogen	Synthesized peptide derived from internal of Human DNL4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total DNL4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA ligase 4; EC 6.5.1.1; DNA ligase IV; Polydeoxyribonucleotide synthase [ATP] 4; LIG4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LIG4

Image



Western blot analysis of extracts from Jurkat cells, using DNL4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.