




TCEB3C Antibody

Product Code	CSB-PA157222
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NG57
Immunogen	Synthesized peptide derived from internal of Human ELOA3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ELOA3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	RNA polymerase II transcription factor SIII subunit A3; Elongin-A3; EloA3; Transcription elongation factor B polypeptide 3C; TCEB3C
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ELOA3

Image



Western blot analysis of extracts from COLO cells, using ELOA3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.