





GADD45G Antibody

Product Code	CSB-PA775807
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95257
Immunogen	Synthesized peptide derived from internal of Human GA45G.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total GA45G protein.
Tested Applications	ELISA, WB, IHC, IF; WB: 1:500-1:3000, IHC: 1:50-1:100, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Growth arrest and DNA-damage-inducible protein GADD45 gamma; Cytokine-responsive protein CR6; DNA-damage-inducible transcript 2; DDIT-2; GADD45G
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GADD45G

Image



Western blot analysis of extracts from HepG2 cells and K562 cells, using GA45G antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using GA45G antibody.



Immunofluorescence analysis of A549 cells, using GA45G antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.