





LAMB1 Antibody

Product Code	CSB-PA264957
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07942
Immunogen	Synthesized peptide derived from C-terminal of Human LAMB1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total LAMB1 protein
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Laminin subunit beta-1 [Precursor]; Laminin B1 chain; LAMB1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LAMB1

Image



Western blot analysis of extracts from HepG2 cells, using LAMB1 antibody.



Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using LAMB1 antibody.



Immunofluorescence analysis of HeLa cells, using LAMB1 antibody.