




MGMT Antibody

Product Code	CSB-PA563873
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16455
Immunogen	Synthesized peptide derived from internal of Human MGMT.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total MGMT protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Methylated-DNA--protein-cysteine methyltransferase; EC 2.1.1.63; 6-O-methylguanine-DNA methyltransferase; MGMT;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MGMT

Image



Western blot analysis of extracts from Jurkat cells, using MGMT antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.