





PHLDA3 Antibody

Product Code	CSB-PA185414
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5J5
Immunogen	Synthesized peptide derived from internal of Human PHLDA3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PHLDA3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Pleckstrin homology-like domain family A member 3; TDAG51/Ipl homolog 1; PHLDA3; TIH1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHLDA3

Image



Western blot analysis of extracts from JK cells, using PHLDA3 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using PHLDA3 antibody.