





TRAF3IP3 Antibody

Product Code	CSB-PA794932
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y228
Immunogen	Synthesized peptide derived from internal of Human T3JAM.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total T3JAM protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TRAF3-interacting JNK-activating modulator; TRAF3-interacting protein 3; TRAF3IP3; T3JAM;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRAF3IP3

Image



Western blot analysis of extracts from 293 cells, using T3JAM antibody.



Immunofluorescence analysis of HeLa cells, using T3JAM antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.