





TP53AIP1 Antibody

Product Code	CSB-PA046843
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCN2
Immunogen	Synthesized peptide derived from internal of Human TPIP1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TPIP1 protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	p53-regulated apoptosis-inducing protein 1; P53aip1; TP53AIP1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TP53AIP1

Image



Immunohistochemistry analysis of paraffin-embedded human thymus gland tissue using TPIP1 antibody.



Immunofluorescence analysis of A549 cells, using TPIP1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.