




POU4F3 Antibody

Product Code	CSB-PA131268
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15319
Immunogen	Synthesized peptide derived from C-terminal of Human POU4F3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total POU4F3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	POU domain class 4 transcription factor 3; Brain-specific homeobox/POU domain protein 3C; Brn-3C; POU4F3; BRN3C
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POU4F3

Image



Western blot analysis of extracts from HuvEc cells, using POU4F3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.