




MRPL14 Antibody

Product Code	CSB-PA962351
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1L8
Immunogen	Synthesized peptide derived from internal of Human MRPL14.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MRPL14 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	39S ribosomal protein L14; mitochondrial; L14mt; MRP-L14; MRP-L32
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MRPL14

Image



Western blot analysis of extracts from HeLa cells, COLO cells and 293 cells, using MRPL14 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.