





RPS12 Antibody

Product Code	CSB-PA831134
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25398
Immunogen	Synthesized peptide derived from C-terminal of Human RPS12.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total RPS12 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	40S ribosomal protein S12; RS12;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RPS12

Image



Western blot analysis of extracts from COLO cells, using RPS12 antibody.



Immunofluorescence analysis of A549 cells, using RPS12 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.