




RPL18 Antibody

Product Code	CSB-PA229371
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07020
Immunogen	Synthesized peptide derived from internal of Human RPL18.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total RPL18 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ribosomal protein L18; RL18;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RPL18

Image



Western blot analysis of extracts from COLO cells and HeLa cells, using RPL18 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.