





RPL34 Antibody

Product Code	CSB-PA299301
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49207
Immunogen	Synthesized peptide derived from internal of Human RPL34.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total RPL34 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ribosomal protein L34; RL34;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RPL34

Image



Western blot analysis of extracts from HepG2 cells, HeLa cells and HUVEC cells, using RPL34 antibody.



Immunofluorescence analysis of HuvEc cells, using RPL34 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.